


## Experiment P-17 Magnetic Field Strength



### Objectives

- To learn about basic properties of magnets.
- To study the relationship between magnetic field strength and the distance from the magnet.

### Modules and Sensors

- PC + NeuLog application
- USB-200 module 
- NUL-214 Magnetic field logger sensor 

### Equipment and Accessories

- |                          |
|--------------------------|
| ▪ Magnet in plastic case |
| ▪ Alnico bar magnet      |
| ▪ Ruler                  |


- The items above are included in the NeuLog Magnets kit, MAG-KIT. The magnet in the plastic case is also included in the NeuLog Mechanics kit, MEC-KIT.

## Introduction


A bar magnet has two ends, known as magnetic poles. One pole is called the north pole of the magnet and the other pole is called the south pole of the magnet. 'Unlike magnet poles' attract each other and 'like magnet poles' repel each other. Magnets also attract other materials which are not normally magnetic (various sorts of metals).

Magnetic fields are force fields produced by electric currents. The magnetic force only act on charged particles in motion. A magnetic field can be represented as lines of force. These lines indicate the direction of the field in a given area, and their number per unit area is proportional to the field strength in that area. The field is stronger near the pole (where the lines converge), and the strength diminishes with the inverse square of the distance from it. The letter ***B*** symbolizes the magnetic field strength and ***l*** symbolizes the distance from the pole. This is an approximate formula in which the results are better as you get closer to the magnet.

$$B \propto \frac{1}{l^2}$$



The Earth has a magnetic field similar to that of a bar magnet. The north pole of a compass needle (a magnet) is attracted to the magnetic south pole of the Earth, which is the geographic north.




In this experiment we will measure the magnetic field strength of a bar magnet as a function of the inverse square of the distance from it.

## Procedure

### Experiment setup


1. Set up the experiment as shown in the picture below.



2. Prepare the two magnets and a ruler.



What happens when you get the two red north poles of the magnets close to each other? What about the two blue south poles? What happens when you get a blue south pole and a red north pole close to each other?

3. Place the ruler in front of the magnetic field sensor. The sensor should be close to the 0 cm line.



4. Place the magnet (with the plastic case) on the ruler. The blue south pole of the magnet should be directed toward the sensor and the magnet's edge should be 14 cm from the sensor.

### **Sensor setup**





5. Connect the USB-200 module  to the PC.
6. Check that the magnetic field sensor  is connected to the USB-200 module.

#### **Note:**

The following application functions are explained in short. It is recommended to practice the NeuLog application functions (as described in the user manual) beforehand.

7. Run the NeuLog application and check that the magnetic field sensor is identified.

## Testing and measurements



8. Click on the **Single Step** icon  to measure the magnetic field strength at 14 cm from the magnet.
9. Click on the **Table** icon  and then on the **Manual values** icon .
10. Insert "Distance" in the column name section. Insert "cm" in the units section. Click on the **Add** button.
11. Insert the relevant distances in the Distance (cm) column after every measurement.
12. Move the magnet closer to the sensor for measuring the strength at 13 cm from the magnet and click on the **Single Step** icon .
13. Continue with the measurements according to the following table.


Measurement number	Distance [cm]
1	14
2	13
3	12
4	11
5	10
6	9
7	8
8	7
9	6
10	5
11	4
12	3

For this measurement, we will assume that the force field lines conversion is at the edge of the magnet.

14. Your data should be similar to the following:



15. Save Click on the **Export** Icon  and then on the **Save value table (.CSV)** button to save your graph.
16. Click on the  icon to go back to the graph.
17. We can see that the magnetic field strength increases as the magnet gets closer to the sensor. The levels increase at a growing rate (exponentially).




18. In order to investigate the relationship between the magnetic field strength and the distance from the magnet, click on the **Functions** icon  and choose Distance as Parameter A.
19. Click on the **Functions** button and choose:  $1/(A^2)$ .
20. A ' $1/(\text{Distance})^2$ ' column should appear on your table.





The screenshot displays the NeuLog software interface in 'Single step' mode. The main window shows a data table with the following columns: Samples, Function 1 (1/(Distance)<sup>2</sup>), Distance (cm), and Magnetic field (mT) ID 1, Exp 1. The data points are as follows:

Samples	Function 1 1 / (Distance) <sup>2</sup>	Distance (cm)	Magnetic field (mT) ID 1, Exp 1
0	0.0051	14	0.141
1	0.0058	13	0.167
2	0.0069	12	0.174
3	0.0083	11	0.188
4	0.01	10	0.228
5	0.0123	9	0.274
6	0.0168	8	0.334
7	0.0204	7	0.413
8	0.0278	6	0.548
9	0.04	6	0.727
10	0.0625	4	1.055
11	0.1111	3	1.588

The interface also features a left sidebar with buttons for 'Function 1' (1 / (Distance)<sup>2</sup>), 'Manual values' (Distance), and 'Magnetic field' (ID 1, Exp 1, mT). At the bottom, there is a toolbar with icons for 'Text size', 'Functions', 'Table', 'Graph', 'Points', and 'Area'.


21. Click on the **Functions** icon  again.
22. Click on the **x axis** button and choose Function 1 as the x axis.
23. Click on the **Graph** icon .
24. Click on the **Zoom fit** icon .
25. Your graph should be similar to the following:



26. We can see that we get a linear relationship between the magnetic field strength and the inverse square of the distance from the magnet.
27. Click on the **Functions** icon  and then on choose Magnetic field as parameter A.
28. Click on the **Function** button and choose Linear fit of A.



29. Click on the **linear fit** module box in the left side of the screen to see the equation.



Function 2  
Linear fit of Magnetic field  
 $13.9904X+0.1063$

## Challenge research

30. Repeat the experiment with the alnico bar magnet.

---

## Summary questions

1. Suppose you have an unmarked magnet. How can you find out where is the north pole and where is the south pole of the magnet?
2. Which magnet is stronger? How stronger is it compared to the weaker magnet?
3. How could you check which magnet is stronger without a sensor?
4. How can you make your own magnet? Point out two ways.