

NeuLog API version 7

This Application Programming Interface (API) specifies how any software should interact with NeuLog sensors.

The API is based on HTTP protocol and can be accessed from any software like C, C++, Python, JS, Java and many more. Also, Microsoft Word™ and Acrobat™ can be used to control the sensors.

Minimum requirements:

- NeuLog sensor(s)
- NeuLog USB module
- NeuLog API software

Getting started:


- Install the NeuLog API software on your PC/MAC/Linux
- Run the NeuLog API
- The API port number will appear on the software screen
- Run your software and contact the NeuLog API via the http "GET" command

The NeuLog API - general:



- The API is based on HTTP protocol and uses the standard “GET” type communication
- Each communication must be initiated by the user
- Each communication looks like this:

[http://localhost:20002/NeuLogAPI?SetSensorRange:\[Light\],\[1\],\[1\]](http://localhost:20002/NeuLogAPI?SetSensorRange:[Light],[1],[1])



- The received command and the result will be shown on the API screen.
- Make sure you are using the right port number.
- The answer is a JSON string such as: `{"SetSensorRange":"True"}`
- This user manual appears when the “help” button ⓘ on the NeuLog API software is clicked.
- Do not use any spaces in the command string.

The API commands:

Command	Description	Parameters	Return
GetServerVersion	Get the server’s software version		<code>{"GetServerVersion":"4.4.4"}</code>
GetSeverStatus	The response could be: Ready, USB missing, Recording		<code>{"GetServerStatus":"Ready"}</code>
GetSensorValue:[sen1.type],[sen1.ID], ..., [seni.type],[seni.ID] Example: GetSensorValue:[Sound],[1],[Light],[2]	Reads one value from the sensor	1. Type & 2. ID of the sensor In the example it is: Sound sensor ID=1 and Light ID=2	<code>{"GetSensorValue":[67.3,345]}</code>
ResetSensor:[sen.type],[sen.ID] Example: ResetSensor:[PH],[2]	Reset the sensor value. Use in force, oxygen and a few other sensors	1. type & 2. ID of the sensor In the example it is: pH sensor ID=2	<code>{"CalibSensor":"True"}</code>
SetPositiveDirection:[sen.type],[sen.ID],[sen.dir] Example: SetPositiveDirection:[Force],[1],[1]	Set the force sensor positive values	1. type & 2. ID of the sensor 3. Direction (1=push positive, 2=pull positive) In the example it is Force sensor ID=1 and the direction is push positive	<code>{"SetPositiveDirection":"True"}</code>

<p>StartExperiment:[sen1.type],[sen1.ID],[sen2.type],[sen2.ID],...[seni.type],[seni.ID],[rate],[samples]</p> <p>Example: StartExperiment:[Temperature],[1],[Temperature],[2],[Temperature],[3],[8],[101]</p>	<p>Start experiment with the specific parameters.</p>	<ol style="list-style-type: none"> Sensors Rate index: <ul style="list-style-type: none"> 1 = 10000 per second 2 = 3000 per second 3 = 2000 per second 4 = 1000 per second 5 = 100 per second 6 = 50 per second 7 = 20 per second 8 = 10 per second 9 = 5 per second 10 = 2 per second 11 = 1 per second 12 = 30 per minute 13 = 15 per minute 14 = 6 per minute 15 = 2 per minute 16 = 1 per minute 17 = 30 per hour 18 = 15 per hour 19 = 6 per hour 20 = 2 per hour 21 = 1 per hour Samples (number) In the example there are 3 temperature sensors; the rate is 10 samples per second, the duration is 101 samples (10 sec) 	<pre>{"StartExperiment": "True"}</pre>
<p>StopExperiment</p>			<pre>{"StopExperiment": "True"}</pre>
<p>GetExperimentSamples:[sen1.type],[sen1.ID],...[seni.type],[seni.ID]</p> <p>Example: GetExperimentSamples:[Light],[1],[Sound],[1]</p>	<p>Get all the samples from specific sensor/s of the current experiment (while running or after running)</p>	<ol style="list-style-type: none"> Sensors 	<pre>{"GetExperimentSamples": [{"Light", 1, 20, 21, 22, 21, 20, ...}, {"Sound", 1, 68.3, 88.8, ...}]}</pre>
<p>SetSensorRange:[sen.type],[sen.ID],[range]</p> <p>Example: SetSensorRange:[Light],[1],[2]</p>		<ol style="list-style-type: none"> Type Id New range 	<pre>{"SetSensorRange": "True"}</pre>
<p>SetRFID:[New ID]</p> <p>Example: SetRFID:[3]</p>		<ol style="list-style-type: none"> New ID 	<pre>{"SetRFID": "True"}</pre>
<p>SetSensorsID:[NewID]</p> <p>Example: SetSensorsID:[3]</p>	<p>Set all connected sensors to the new ID</p>	<ol style="list-style-type: none"> New ID <p>In the example the new ID is 3</p>	<pre>{"SetSensorsID": "True"}</pre>

Few examples working with Sound and Light sensor:

<http://localhost:22001/NeuLogAPI?GetServerVersion>

<http://localhost:22001/NeuLogAPI?GetSeverStatus>

[http://localhost:22001/NeuLogAPI?GetSensorValue:\[Sound\],\[1\],\[Light\],\[1\]](http://localhost:22001/NeuLogAPI?GetSensorValue:[Sound],[1],[Light],[1])

[http://localhost:22001/NeuLogAPI?StartExperiment:\[Sound\],\[1\],\[Light\],\[1\],\[8\],\[101\]](http://localhost:22001/NeuLogAPI?StartExperiment:[Sound],[1],[Light],[1],[8],[101])

<http://localhost:22001/NeuLogAPI?StopExperiment>

<http://localhost:22001/NeuLogAPI?GetExperimentSamples>

[http://localhost:22001/NeuLogAPI?SetSensorRange:\[Light\],\[1\],\[2\]](http://localhost:22001/NeuLogAPI?SetSensorRange:[Light],[1],[2])

Sensor list:

'Temperature', 'Light', 'Voltage', 'Current', 'PH', 'Oxygen', 'PhotoGate', 'Pulse', 'Force', 'Sound',
'Humidity', 'Pressure', 'Motion', 'Magtnetic', 'Conductivity', 'GSR', 'CO2', 'Barometer', 'Rotary',
'Acceleration', 'Spirometer', 'SoilMoisture', 'Turbidity', 'UVB', 'EKG', 'Colorimeter', 'DropCounter',
'FlowRate', 'ForcePlate', 'BloodPressure', 'Salinity', 'UVA', 'SurfaceTemp', 'WideRangeTemp',
'InfraredThermometer', 'Respiration', 'HandDynamometer', 'Calcium', 'Chloride', 'Ammonium', 'Nitrate',
'Anemometer', 'GPS', 'Gyroscope', 'DewPoint', 'Charge'