

## Experiment S-1

# Light in Different Environments

## Objectives






- To learn about qualitative and quantitative data
- To explore Light changes in different environments

## Description


In this experiment you will measure light intensity in different environments, indoors and outdoors, through Panda's light sensor.



## Instructions

- Press on the icon on the top left of the Panda screen and pick the "Light" option (if it is not already selected).
- Press on the  icon and change the range to 0-6000 lx.
- Measure light in the following conditions:
  - Indoors with no light and the windows closed
  - Indoors with no light and the windows open
  - Indoors with the windows open and lights on.
- Write the light intensity in the next page.
- Press on the  icon and change the range to 0-240,000 lx.
- Measure light in the following conditions:
  - Outdoors (places the sensor horizontally)
  - Outdoors with the sensor aimed at the sun
- Write the light intensity in the in the next page.
- Observe the Light change in other Panda modes: Column , Graph , Gauge 





## Light Intensity Indoors


Express the light intensity in words (Qualitative data):

Measured value (Quantitative data): :

No light and windows closed

No light and windows open

With the windows open and lights on



## Light Intensity Outdoors

Express the light intensity in words (Qualitative data):

Measured value (Quantitative data):

Sensor placed horizontally

Sensor aimed at the sun

What are the advantages in quantitative and qualitative measurements?

---



---



---



---



---



---

