

# Sense of Touch




## Experiment



### Objectives

- Learn about the **skin** and how objects touching it affect sensory perception.
- Identify which parts of the body are the **most sensitive to touch**.
- To understand that the **skin is the largest sensory organ in the body**.
- Explore how the skin detects touch, pressure, temperature, and pain.
- Examine the **three major layers of the skin**.
- To recognize that **fingers, lips, genitals, and toes** are among the most **sensitive areas** of the body.

### Modules and Sensors

- PC + NeuLog application
- USB-200 module  or BLT-202 Bluetooth module 
- NUL-217 GSR (Galvanic Skin Response) logger sensor 

### Equipment and Accessories

▪ Wash bottle (or a glass)
▪ Cotton pad
▪ Pen

### Materials

▪ Water for wash bottle	–
-------------------------	---

\* This experiment is based on Morris Tischler's book 'Concepts of Biotechnology'.

## Introduction

Emotional and sensory stimulation **trigger** the production of sweat. The GSR (Galvanic Skin Response) logger sensor measures sweat secretion from the sweat glands in the hand. As stimulation occurs, the amount of sweat **increases**, leading to **higher skin conductivity**.

### Why does emotional and sensory stimulation cause sweating?

The **fight-or-flight response** is activated in reaction to an **acute survival threat**, preparing the body to **react or retreat**.

When a threat is perceived, the **sympathetic nervous system** (a part of the autonomic nervous system) is activated, leading to the **release of noradrenaline and adrenaline**. These hormones bind to **adrenergic receptors** in peripheral tissues, triggering:

- **Pupil dilation**
- **Increased heart rate**
- **Elevated blood pressure**
- **Faster breathing**
- **Sweat production**

Sweat is released during the fight-or-flight response to **help dissipate excess heat** generated by increased muscle activity.

Humans have **several millions of sweat glands**, located in the **middle layer of the skin (dermis)**. Sweat is transported to the **surface layer** (epidermis) through **sweat ducts**. Sweat glands are concentrated in areas such as the forehead, palms, , armpits, and soles of the feet.

## THE SKIN AS a Sensory Organ


The skin is our largest sensory organ in the human body, **covering** its entire **surface**. In a person standing 1.8 meters (about 6 ft.) tall, the skin measures approximately 1.75 sq. meters (3,000 sq. inches) and accounts for about 15% of total body weight.

The skin is not only a sensory organ but also serves as a protective **barrier** that insulates the **body's** nerves, muscles, and organs from mechanical and electrical injury. Additionally, it helps retain vital body fluids. The skin protects the body from the sun's harmful **ultraviolet** (UV) and infrared radiation and regulates temperature to prevent excessive fluctuations.

The skin contains a variety of **nerve endings** that are **highly responsive** to pressure, pain, and **temperature changes (hot or cold)**.

**Continuous pressure applied to the skin, particularly on the fingers, can lead to a temporary reduction in sensitivity.**

The following figure illustrates the **most and least sensitive** of the body. **The skin** around the mouth, fingers, and toes is the most sensitive, as these areas contain the **highest concentration** of nerve fibers.



The skin is sensitive to several **types** of stimulation:

1. Touch
2. Pressure
3. Cold
4. Heat
5. Pain

The fingers **contain up to** 800 nerve endings per square centimeter (1 cm x 1 cm). The following figure illustrates the areas of greatest sensitivity in the hand.



**Sensitive Areas of the Hand**

#### **GSR Measurement ranges:**

The NUL-217 NeuLog GSR logger sensor has two measurement modes:  $\mu$ S (**microSiemens**) and Arb (**arbitrary units**).

- The  **$\mu$ S range** measures **absolute conductivity** in microSiemens, which varies between individuals.
- The **Arb range** amplifies changes in conductivity, making it more useful for detecting **relative responses**.

When a person's **baseline conductivity is high**, **small changes may not be easily detected**.

Therefore, **Arb mode is preferred when tracking changes is more important than measuring absolute values**, such as in **lie detection experiments**.


## Procedure

### Experiment setup

1. Set up the experiment as shown in the picture below.





2. Ensure you have the following materials:
  - a wash bottle filled with water (or a glass of water),
  - a cotton pad,
  - five blank cards, a pen
  - a bottle of alcohol, perfume or cologne.
3. Turn off the air-conditioning system in the room to minimize external influences on the results.
4. **Note:** You cannot measure **your own** emotional responses.
  - Two students will conduct the experiment.
  - A third student will be the test subject.
5. Instruct the subject to **warm their hands** by rubbing them together.
6. Have the subject **moisten** two fingers using a wet cotton pad.
7. Attach the **finger electrodes** to the moistened areas as shown in the figure.



8. Instruct the subject to sit down.  
**Stand or sit behind them** to ensure they **cannot see you, the other student, or the computer screen.**

## Sensor setup

9. Connect the USB-200 module  to the PC.
10. Ensure that the GSR sensor  is properly connected to the USB-200 module.


### Note:



The following application functions are **briefly explained**. It is recommended to practice using the NeuLog application in advance by following the **user manual**.

11. **Launch** the NeuLog application and verify that the GSR sensor is detected.

## Settings

11. Click on the **Sensor's Module** box.
12. Select the **Arb mode** to change the sensor's measurement setting. .
  - The **μS (microSiemens) mode** measures **absolute conductivity**, which varies between individuals.
  - The **Arb (arbitrary units) mode** amplifies changes in conductivity, making it better for detecting **relative responses**.





13. Click on the  icon to return to the graph.
14. Click on the **Run Experiment** icon  and set the following parameters:
  - **Experiment duration:** 10 minutes
  - **Sampling rate:** 10 samples per second

## Stimulating response

People **rely heavily** on their **sight and hearing**, so **minimizing these senses** enhances attention to **other sensory inputs**. Touch and smell **stimulation** is more effective **without visual or auditory distractions**, as these sensory become more sensitive in their activity.

15. Instruct the subject to **sit quietly** for a few moments with eyes closed, head down, legs uncrossed, and arms resting on the armchair supports or lap.
  - **Ensure the subject does not exert pressure** on the finger electrodes.
16. **Maintain silence** to prevent auditory stimulation. These conditions should be maintained throughout the measurement.
17. Observe the measured response value in **Arb units** on the sensor.

**Note:** This value varies between subjects and is influenced by **emotional state**.

18. Click on the GSR sensor module box and set the **Y-axis limits** to be approximately 6000 **units** above and below the measured value.
19. One student will **perform the stimulations**, while the other will **record the time of each event** and the **corresponding Arb** value.
20. Click on the **Record** icon  to start the measurement.
21. Click on the **Arrows** icon  to view the sensor's values during the measurement.
22. The graph may **fluctuate**. Wait silently until only **small variations** appear.

## **Touch Sensitivity Testing**

**Event 1:** Lightly touch the subject's ear and observe the response on the graph.

- The production of **even** a small amount of sweat will **alter** the measured value.
- **More sweat production** → **Greater change in conductivity**.

**Note:**

**There is a time delay (0.8 to 4 seconds) etween stimulation and the sensor's response** as sweat is produced.

**Event 2:** Lightly touch the other ear and observe the response.

**Event 3:** Blow gently across the subject's ear or face after repeated touch stimulation to see if they respond to air movement.

- **Habituation: A decrease in response** due to repeated exposure to the same stimulus.

## Event 4: Conditioning Test

Conditioning is the process of training a subject to perform in a certain way.

**Touch the subject's skin** in each of the following areas and **observe the graph**:

- Shoulder**
- Calf** (lower leg)
- Forearm**
- One toe**
- One toe** → Wait → **Then touch the arch of the foot**
- Forefinger (index finger)**
- Forehead**
- Nose, then lips**

23. **Repeat the experiment** with another subject, and compare the results. Observe how **different individuals respond differently** to touch stimuli.

### Note:


- This experiment does **not exclusively measure touch sensitivity**.
- **Emotional factors** can also **influence the results**.
- Simply **thinking about a frightening event** may **trigger a response** on the recording device.

## PAIN MEASUREMENT

24. Observe how the subject responds to **mild pain stimuli**.
  - Use a **toothpick or pin** to touch the areas that were previously tested.
  - Apply **light pressure** to induce a **very slight pain sensation** at the contact point.
  - **Compare the response** to previous measurements taken **without pain stimulation**.
  - **Is the response greater than that recorded for non-painful touch?**

### NOTE:

- Touch the skin quickly rather than pressing slowly.
- **Why?** Slow, sustained pressure allows the skin to **adapt to pain**, reducing the response.

25. Click on the **Export** Icon  then select **Save value table (.CSV)** to save your graph.



## Conclusions

1. Emotional and sensory stimulation (through the five classic senses) triggers sweat production, which is a key component of the **fight-or-flight response**.
2. The graph confirms that when the subject is **under stress**, sweat secretion **increases** and can be accurately measured using the GSR logger sensor.
3. Certain areas of the body have **higher concentrations of sweat glands**, such as the forehead, palms, armpits, and the soles of the feet. Sweat measurements in regions like the forehead, back or chest are more **challenging** to obtain due to **variability in gland distribution**.
4. Sweat level **changes** occur **within** 0.8 to 4 seconds after the subject experiences stimulation.


## Summary questions

1. Based on your results, which **areas** of the body **exhibited the highest sensitivity to touch**?

## Appendix

### SKIN STRUCTURE

The skin consists of multiple layers of cells, with sensory nerve endings located beneath the outer layers. The following figure provides a simplified illustration of the skin's structure



#### Sensory Receptors in the Layers of the Skin

The **three major layers of the skin** are:

1. **Epidermis (outer layer)**
2. **Dermis (inner layer)**
3. **Subcutaneous layer (deepest layer, containing fatty tissue)**

Hair fibers which are anchored in the subcutaneous **layer**, penetrate through the epidermis.

- Near the root of each hair, nerve endings around the follicle.
- If the hair is pulled, these nerve endings send a **pain signal** to the brain.

The epidermis **undergoes continuous renewal**:

- **Cells migrate** to the surface over **1–3 months**, replacing older cells that **shed naturally**.
- **Melanocytes** in the **deepest epidermal layer** produce **melanin**, which determines **skin color**.

#### Factors Affecting Skin Color

- Melanin quantity, size, and distribution influence skin tone, not race or sex.
- Sun exposure triggers melanin production, resulting in tanning.
- Carotene, a yellow pigment, also contributes to skin coloration.
- Albinos, who lack melanin production, have very pale skin and are highly sensitive to UV radiation.

## Touch & Sensory Receptors in the Skin

The skin contains **various sensory receptors**, each specialized for detecting different stimuli:

1. **Free nerve endings** – Found throughout the skin and in sensory organs like the eyes, responsible for **light touch and pressure sensitivity**.
2. **Meissner's corpuscles** – **Encapsulated nerve endings** that are **highly sensitive to touch and low-frequency vibrations**. These are abundant in the **fingertips, lips, and genitals** and adapt **within seconds**.
3. **Merkel's discs** – These receptors **maintain a continuous response to steady pressure**, allowing the brain to **detect sustained contact** (e.g., clothing against the skin).
4. **Hair follicle receptors** – Nerve fibers **wrapped around hair roots**, responding to **even slight hair movement**.
5. The **brain becomes aware of the burn** only **a few seconds later**.
6. This rapid **reflex arc mechanism** is **essential for survival**, preventing **serious injuries** before the brain even processes the pain.

Another type of **sensory receptor** is the **Pacinian corpuscle** which detects **rapid changes in pressure** and **vibrations**.

- When pressure is applied, **fluid shifts within the corpuscle**, generating **vibrations** that are transmitted to the **brain** as touch and pressure sensations.

The **immune system** plays a vital role in **skin protection**:

- It sends signals via the lymphatic system to **detect** foreign **substances** in contact with the skin.
- It reacts rapidly to **objects piercing the skin or causing injury**, triggering **inflammatory and immune responses**.

As shown in the diagram, different nerve fibers are distributed at **various depths** beneath the skin.

- The nerve fibers **respond to** surface changes, making some areas **more sensitive** than others.

Highly sensitive areas include:

- Fingertips
- Mouth region
- Tip of the tongue
- The triangular area between the mouth and lower nose is particularly sensitive to touch.

- Less sensitive areas include:

- legs,
- arms,
- trunk of the body.

The brain contains 99% of the nervous system, yet the remaining 1% controls most involuntary body functions without conscious thought.

## THE REFLEX ARC

The nervous system consists of two main pathways:

1. Voluntary actions - Controlled by the brain and its subsections.
2. Involuntarily actions - Controlled by the spinal cord, bypassing conscious thought.

Some actions require an immediate response without brain involvement.

### Example: The Knee-Jerk Reflex

- Tapping the knee with a **rubber mallet** causes an **automatic leg movement**.
- This involuntary response is called a reflex arc, defined as a rapid **automatic reaction** to a stimulus.
- The autonomic nervous system also controls many routine functions, including:
  - heartbeat **regulation**,
  - breathing patterns,
  - Skin temperature control.

The following figure **illustrates** a person **reflexively pulling** away from a painful stimulus without **conscious thought**. This is an example of how the **reflex arc** functions.

1. The skin receptors **detect the painful stimulus (e.g., heat)**.
2. A rapid **nerve signal** (action potential) is transmitted through the **afferent nerves (sensory neurons)** to the spinal cord.
  - Afferent nerves carry signals **towards** the spinal cord or brain.
3. The **spinal cord processes the information and immediately sends a response** via the **efferent nerves (motor neurons)**.
  - **Efferent nerves** carry signals **away** from the spinal cord or brain to the muscles.
4. The **effector muscle contracts**, causing the person to **quickly withdraw from the stimulus**—all within **milliseconds**.

This **reflexive response** occurs **before the brain fully processes the pain**, ensuring **rapid protection from harm**.

